

VIEW THE RADIOLOGY READING STATION FOR VETERINARIANS

VET READINGS MADE SIMPLE



VERSION 2.7

VET-VIEW



VET-VIEW 2.7 SPECIALIZED FOR VETERINARY NEEDS

The most important issue for vets who frequently read imaging studies is selecting the appropriate reading station.

VET-VIEW is an easy-to-use reading station, specialized for veterinarians. A customized veterinary graphical user interface and the especially adapted text overlay are only some of the many veterinary features. All important veterinary data, like breed description or studbook name, is visible on every image, thus providing easy patient recognition. VET-VIEW is available in two different versions: BASIC and PRO. The BASIC version includes all required features for viewing and reading while the PRO version offers professional tools for reading specialists like presentation states and individually customizable hanging protocols. Hanging Protocols make it easy to view for instance pre-purchase examinations of horses for facilities dealing with equine patients. Studies and images can be loaded automatically as each veterinarian wishes to read them.

Vet-specific measurements, such as the Norberg angle, make image review simple and fast.

Due to its special software design, VET-VIEW is very fast. An automa-

VET-VIEW

tic tiling function arranges imaging studies on the screen without any user interaction. A configurable toolbar in the viewer allows the easy selection of tools that are needed for image processing. The most relevant buttons can simply be set up in the personal default toolbar.

VET-VIEW can be integrated into virtually any practice management system, EMR or HIS. Customized modules for 3D post-processing, nuclear medicine, orthopedic templating, image stitching and more equip specialty readers with the required tools.

VET-VIEW even contains an integrated basic structuring reporting module for creating, editing, storing and sending individualized structured reports.

It is further possible to write imaging data and reports on CDs, DVDs or USB media, to print on Windows® laser printers or DICOM medical imagers and to email any study to a colleague - even in diagnostic quality.

VET-VIEW may be used as a radiological media burning station – or rather, as a semi-automatic burner/ printer in connection with iQ-ROBOT. VET-VIEW can also work as an image acquisition station by importing images from scanners, directories and portable media, including options for patient reconciliation offered in the PRO version. Scanned documents can simply be attached to a study.

A special feature of the import module are the fields for veterinary view positions, body parts and vetspecific information, such as species, breed, owner and studbook info, which can be added to the image to complete the information workflow. VET-VIEW is ideal for teleradiology purposes. Encrypted DICOM email makes it possible: Ask any colleague anywhere in the world for a second opinion by just doing two mouseclicks! This new technology offers easy and fast information exchange among colleagues by email. Or use the email feature to send out image data as JPEG files.

The application is compatible with virtually all medical imaging DICOM devices on the market, regardless of the manufacturer.

THE BEST FOR VETS



VET-VIEW SCREENSHOTS




According to veterinary DICOM standard:				
owner's last name and animal name in one search fiel				
Fallaci = owner's last name Capricciosa = animal name				
		/		
Animal ID	Owner / Animal name	Referring veterinarian		
	Fallaci, Capricciosa 🗡			
Accession no.	Date of birth	Study description		
_				
From: 14.03.2012	To: 14.03.2012	Status		



- A The Study Browser provides a clear overview of all imaging data including structured reports. A patient root view will automatically sort all available studies by patient.
- All image processing tools can be accessed either by selecting buttons or menu items or even by using shortcuts. For an easy orientation within a series, the userfriendly viewer offers several scout pilot tools, such as the lines mode and the 3D position display.
- **C** By using the search filters it is possible to narrow down the search for intended studies. This can save a lot of time when querying a PACS with a high volume of stored data.
- VET-VIEW can be used as a printing station; the integrated print manager allows for a flexible optimization of the printouts.








- E Single and multiple studies and series can be exported to a CD, DVD or memory stick, all including a VET-LITE CD viewer. Exporting to DICOM email, to directories and different file formats is also possible.
- The import module is a comprehensive tool for the acquisition of images from different sources (such as scanners, capture devices, cameras, microscopes or file directories) and their conversion to DICOM.

VET-VIEW WORKFLOW



VET-VIEW places the reading veterinarian at the heart of all diagnostic processes: Studies from all modalities or a PACS are received in the reading station, colleagues and referring veterinarians can be connected via DICOM or DICOM email and VET-VIEW can easily be integrated with any practice management system (PMS), EMR or HIS. Even specialty tools for 3D post-processing, orthopedic templating, nuclear medicine or stitching are closely interfaced with the viewer.

VET-VIEW FEATURES

GENERAL	• 64 bit support for Windows 7 Professional or bigher	
GENERAL		
	 Scalable buttons and possibility to set font size (for high-resolution displays) 	
	 Multimonitor support 	
STUDY BROWSER	Study table	
	 Thumbnail preview 	
	 Virtual patient root model 	
VIEWER	 Veterinary text overlay 	
	Image processing (including comparison of multiple studies, image filters and	
	synchronization/desynchronization of series)	
	 Configurable bottom toolbar and creation of viewer shortcuts 	
	 Measurements (including distance, angle, Cobb angle, Norberg angle, ROI, 	
	modify/delete options)	
	 Basic annotation of images (incl. Secondary Capture image creation) 	
	 DICOM presentation states on image level 	
	 Configurable scoutlines to display plane projections and/or plane intersections 	
	Integration of 3rd party image processing or image analysis tools	
	 Support of DICOM shutter sequences 	
REPORTING	Structured reporting module with read/write/edit support of Basic SR objects	
	 Read-only support for other SR formats 	
IMPORT	Lists containing veterinary body parts and image views	
	 Additional veterinary DICOM information included, e.g. studbook, breed, species or 	
	owner information, sex neutered and more	
	 TWAIN interface (scanner, video capture and CR) 	
	 Drop box interface for automatic image file import 	
	 Image file import (BMP, JPEG, TIFF, RAW) 	
	 Patient CD/DVD/USB import 	
EXPORT	- Batiant CD/DVD/USP	
	 Patient CD/DVD/USB Integration of direct interface to iQ-ROBOT module for patient CD/DVD creation and 	
	label printing (not available in the USA)	
	 Image export to image file or AVI video file (including measurements, windowing, 	
	color remapping, annotations and more)	
	 Sending of JPEG images by email 	
	- Sending of Stield images by email	
PRINTING	Print to DICOM imagers of any brand	
	 Windows[®] print 	
	 True-size (1:1) DICOM printing 	
COMMUNICATION	Server administration tool for easy configuration	
	 Smooth interoperability with iQ-products; support of various transfer syntaxes 	
	 HIS/PMS interface (accession number or GDT/BDT) (others on request) 	
	 DICOM Query/Retrieve, C-FIND, C-STORE, DICOM print client and DICOM email 	

VET-VIEW FEATURES

ADD-ONS	 Optional 3D post-processing add-on module (iQ-3D) Optional image stitching software (iQ-STITCH) Optional orthopedic templating module (OrthoView[™]) Optional nuclear image processing (iQ-NUC) 	
LICENSING	 Single (PC-based) or concurrent (floating user) licenses 	
LANGUAGES	 English, French, German 	

VET-VIEW PRO FEATURES

IN ADDITION TO VET-VIEW FEATURES

VIEWER	 Hanging Protocols (including the creation of HP sequences and modality-specific toolbars) Perpendicular distance, polygonal ROI Copy of measurements to other images to enable 3D measurements Dynamic magnifier window for seamless image zooming Annotation of images (incl. storing and sending as Presentation States) Color schemes for nuclear medicine, e.g. GE Color, Hot Iron, Lut7
COMMUNICATION	 DICOM Worklist Client
	 Video capture interface (DirectShow®)* Foot switch support

• VET-CR ACE post-processing incl. histogram function and S-value (optional)



* We recommend the VET-CAPTURE or VET-CAPTURE PRO hardware packages for an efficient and guaranteed use. They each include a frame grabber card and a foot switch.

HARDWARE & SOFTWARE REQUIREMENTS				
	мілімим	RECOMMENDED		
OS:	Windows XP Pro 32 bit Windows 7 Pro 32 bit or 64 bit	Windows XP Pro 32 bit Windows 7 Pro 32 bit		
CPU:	Intel [®] CPU with 1 GHz min. 1.5 GHz for cardiology	Intel® Dual-Core CPU > 2 GHz		
RAM:	1 GB RAM min. 2 GB RAM for cardiology	4 GB RAM		
HDD:	2 GB free hard disk space (DMA33)	20-500 GB free hard disk space (S-ATA-II)		
Network:	10 Mbit/s network adapter	100 Mbit/s network adapter		
Graphics:	24 bit color or 8 bit gray output Any NVIDIA, ATI or Matrox graphics card	24 bit color or true 12 bit gray output iQ-GRAFIX or iQ-GRAFIX PRO		
Display:	1024 x 768 pixel High-resolution monitors for diagnostics	Dual screen of 1280 x 1024 pixel or more IMAGE DISPLAYS diagnostic monitors		
Peripherals:	CD-ROM drive	Scroll mouse, CD or DVD writer, PostScript printer		
Hardware:		Dell		

OUR PRODUCTS FOR YOUR IMAGING NEEDS

VET-PACS VET-VIEW VET-RIS IMAGE DISPLAYS The full-featured, reliable and affordable PACS for vets The radiology reading station designed for veterinarians The radiology information system for veterinarians Medical diagnostic monitors

IMAGE INFORMATION SYSTEMS LTD. OFFICE 404, 4TH FLOOR | ALBANY HOUSE | 324/326 REGENT STREET | LONDON W1B 3HH | UNITED KINGDOM TEL. UK: +44 207 193 06 20 | TEL. GER: +49 381 496 582 0 | TEL. US: +1 704 323 66 63 FAX UK: +44 207 976 48 97 | FAX GER: +49 381 496 582 99 | FAX US: +1 704 945 71 01 WWW.IMAGE-SYSTEMS.BIZ | INFO@IMAGE-SYSTEMS.BIZ